
Tissue-Tek®
Block Filing Cabinet 2400

*Reliable storage with
superior capacity*



continuous innovation for pathology



Modern ergonomic design providing best-in-class 2,400 block capacity for convenient and easy storage

Tissue-Tek®
Block Filing Cabinet 2400



continuous innovation for pathology



Tissue-Tek®
Block Filing Cabinet 2400



continuous innovation for pathology



Tissue-Tek®
Block Filing Cabinet 2400



The Tissue-Tek® Block Filing Cabinet 2400 provides an efficient solution for organising and storing paraffin-embedded specimen blocks, featuring a high-capacity and ergonomic design.

- Each cabinet stores 2,400 individual blocks
- Includes adhesive and removable drawer labels for easy information recording and reusability
- Each label features a unique barcode ID for efficient electronic tracking
- Deep drawers with angled label surface ensure easy block placement and lookup while minimising user strain
- Ergonomic exterior features include a recessed base, while interlocking surfaces allow for safe moving and stacking of cabinets
- Stackable up to three units, for high-density storage
- Latch mechanism prevents drawers opening accidentally during transport
- High-grade and durable plastic prevents storage degradation in humid or harsh environmental conditions
- Includes foam inserts to keep blocks upright during filling

Specifications



Product name	Tissue-Tek® Block Filing Cabinet 2400
Item code	4193-E
Applications	Storage of formalin-fixed paraffin-embedded specimen blocks
Number of drawers	6
Capacity	2,400 paraffin-embedded specimen blocks
Dimensions	43.7 (W) x 27.0 (D) x 42.0 (H) cm
Weight	Empty: 6.8 kg Filled: ~21.3 kg
Stackable	up to 3 units

Accessory list

Amount in packaging	3
Item number	4197
Item description	Tissue-Tek® Block Filing Cabinet Drawer Label Set



Visit us at sakura.com